www.curaden.ch f4r du veta mer om den nya schweiziska munv4rddsskolan.

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**GB How to clean my teeth properly**

Join the first animated, practical instruction according to the "New Swiss School of Oral hygiene" for a correct use of the toothbrush, floss and interproximal brushes. **Visit the [www.curaden.ch](http://www.curaden.ch) homepage to learn how to maintain fresh breath, healthy gums and teeth without plaque and stain.**

1. Hold your CURAPROX Single Sulcular brush at a right angle to your tooth, inserting the outer edge of the brush carefully against the gum line.
2. Hold your single sulcular brush at an angle against the gum.
3. Insert the outer edge of the brush very carefully against the gum line.
4. Slowly guide your single sulcular brush along the gum line with very fine vertical movements.
5. The brush must be kept in a vertical position against the tooth you are cleaning at all times.
6. The cleaning action should start on one side of your tooth and finish on the opposite side of the same tooth.
7. Continue cleaning the entire lingual surface.
8. Carefully feel your way along the gum line and clean your tooth according to the front-side.
9. The cleaning action will end on the opposite side of your tooth.

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